



MN-4

Coarse Manipulator

Manual coarse manipulator with many convenient functions.

This manual coarse manipulator is used in combination with a hydraulic micromanipulator to transfer the pipette into the microscope field of view. Since its operation is manual, this model offers excellent response at affordable price. Compared with the MMN-1, setup and installation are straightforward. Rotating mechanism is included for easy specimen or pipette exchange.



Specification

Accessory included		Hex Wrench
Movement range	Coarse	X30mm, Y30mm, Z30mm
Dimensions/Weight		W125 × D110 × H135mm, 600g