



# IM-HD1

## Pneumatic Injector (for Holding Side)

### Injector for the holding side attachable to the manipulator joystick

This pneumatic injector is intended exclusively for holding an egg. It can be attached to the manipulator joystick, which allows holding and releasing an egg in a delicate and secure manner operated only with a forefinger. It enables doing micromanipulation and suctioning work at the same time. Two types of tubing which has different inner diameter are included to customize the response speed.



\* IM-HD1 can be attached to the joystick of [MMO-4](#), [MMO-202ND](#), [MO-202U](#), [MOM-202D](#), [MON-202D](#), [MWO-202D](#), [ON](#) series and [NT](#) series.

\* When attaching to [TAKANOME™](#), refer to [IM-HD1T](#).

## Specification

<b>Accessories included</b>	Injection Holder ( <a href="#">HI-9</a> ) Cleaning Rod Tubing ( <a href="#">CT-4</a> ) Tubing (CT-1-1.5) Silicone Rubber Gasket (HI01PK01) Hex Wrench
<b>Dimensions/Weight</b>	W24 × D55 × H54mm, 58g

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2018