

# ***NARISHIGE***

## **G-1/1.2/1.5/2/3** **Glass Capillary**

### **General type of glass capillary.**

ø1mm ~ ø3mm O.D. glass capillaries are available for forming various shapes of electrodes. G-1, G-1.2, and G-1.5 are pre-washed in an ultrasonic washing machine.



# NARISHIGE

## G-1/1.2/1.5/2/3 Glass Capillary

### Specification

Type	G-1	G-1.2	G-1.5	G-2	G-3	GD-1	GD-1.2	GD-1.5	GD-2
Outer diameter (mm)	1	1.2	1.5	2	3	1	1.2	1.5	2
Inner diameter (approx. mm)	0.6	0.75	0.9	0.8	1.8	0.6	0.75	0.9	1.2
Length (mm)	90*	90*	90*	90*	90	90*	90*	90*	90*
Contents (pcs)	500	500	500	200	200	500	500	500	200
Material	Borosilicate								

Type	GC-1	GC-1.2	GC-1.5	GDC-1	GDC-1.2	GDC-1.5	G-100	G-100L	G-1000
Outer diameter (mm)	1	1.2	1.5	1	1.2	1.5	1	1	1
Inner diameter (approx. mm)	0.6	0.75	0.9	0.6	0.75	0.9	0.75	0.75	-
Length (mm)	90	90	90	90	90	90	90*	100*	90*
Contents (pcs)	200	200	200	200	200	200	500	500	200
Material	Borosilicate								

\* A custom length between 70mm and 180mm is also available by special order.

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2020