



IM-HD1T

Pneumatic Injector (for Holding Side)

Injector for the holding side exclusively for TAKANOME™

This pneumatic injector is intended exclusively for holding an egg. It can be attached to the joystick of TAKANOME™, which allows holding and releasing an egg in a delicate and secure manner operated only with a forefinger. It enables doing micromanipulation and suctioning work at the same time. Two types of tubing which has different inner diameter are included to customize the response speed.



* IM-HD1T can be attached to [MTK-1](#) series.

* When attaching to the joystick of conventional models, refer to [IM-HD1](#).

Specification

Accessories included	Injection Holder (HI-9) Cleaning Rod Tubing (CT-4) Tubing (CT-1-1.5) Silicone Rubber Gasket (HI01PK01) Hex Wrench
Dimensions/Weight	W24 × D55 × H54mm, 58g

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2018